

FEATURES

MM5LH

- ◆ Digital LCD Display.
- ◆ Last Measurement in Memory.
- ◆ Heating & Stirring can proceed Simultaneously.
- ◆ Bushless Motor for Low Noise.
- ◆ Anti scratch and Corrosion Resistance.
- ◆ Ceramic Plate Surface.
- ◆ ABS Cabinet.



SPECIFICATIONS

HSN : 84798200

PRODUCT NAME	CERAMIC HOTPLATE MAGNETIC STIRRER
MODEL	MM5LH
Speed range	100-2000rpm
Maximum mixing capacity	5000 ml
Temp. Range & Accuracy	300°C Accuracy $\pm 1^{\circ}\text{C}$
Temperature Sensor	PT100
Time range	0-99h.99min
Top Plat Size	168 x 168mm
Display	LCD
Motor Type	Bushless Motor
Operating Mode	Temperature, Speed, Timing, Memory (last data)
Power	500W
Power Supply	AC 220V / 50-60Hz
Packing Dimension (LxWxH)	335 x 235 x 170mm
Net / Gross Weight	2.5kg / 3kg

* Specifications are subject to change without Notice