

FEATURES **MWB4, MWB6**



MWB4



MWB6

- ◆ Digital Thermostatic Display.
- ◆ Stainless Steel Tank.
- ◆ High Precise Temperature Control.
- ◆ High Corrosion Resistance.

SPECIFICATIONS

HSN : 84192090

PRODUCT NAME	DIGITAL WATER BATH	
MODEL	MWB4	MWB6
Column	2 Rows 4 Holes	2 Rows 6 Holes
Internal volume (approx.)	12L	18L
Temp. Range	RT+5~99.9°C	
Temp. Resolution, Accuracy	0.1°C, Accuracy ±0.5°C	
Water sink	Stainless steel	
Temp. control mode	PID	
Cover aperture	120mm	
Cover Ring per hole	5pcs.	
Power Supply	AC 220V/ 3.6A(50/60Hz)	
Power rating	800W	1200W
Inner Chamber size (LxWxH)	300 x 300 x 140 mm	450 x 300 x 130 mm
Packing Dimension (LxWxH)	410 x 410 x 260 mm	660 x 400 x 260 mm
Gross Weight	6kg	8kg

* Specifications are subject to change without Notice