

FEATURES

MIM20C

- ◆ Ice Crushing Facility.
- ◆ Housing Material Stainless Steel Type.
- ◆ R290A Refrigerant.
- ◆ Water Inlet Manual & Automatic.
- ◆ Bullet Shape Ice.
- ◆ Easy to move through hand grip.
- ◆ Low Water Alarm.
- ◆ Ideal Working Temperature 20°C.
- ◆ Transparent Window on Top Cover to view inside.



SPECIFICATIONS

HSN : 84186990

PRODUCT NAME	ICE MAKER & CRUSHERS
MODEL	MIM20C
Ice capacity	20kg/24H
Function	Ice making / Crushing ice
Storage Bin	3kg
Refrigerant	R290A
Ice Shape	Bullet Shape
Housing material	Stainless Steel type
Water inlet method	Manual / Automatic Water inlet
Ice Output Each Cycle	18 pieces
Ice Making Time	10-12 minutes in Ideal Conditions
Input power	180W
Power supply	AC 220V/50Hz
Packing Dimension (L x W x H)	445 x 420 x 465mm
Net / Gross Weight	15.6kg / 17.2Kg

* Specifications are subject to change without Notice

Monotech Lab Instruments and Calibration

B-9, Shabnam Building, MEK Industrial Estate, Opp. Fire Brigade, Mumbra, Dist. Thane 400612 Maharashtra