

FEATURES

MIF30D, MIF50D, MIF100D

- ◆ RED LED Display.
- ◆ Fully Automatic microprocessor control.
- ◆ Stainless Steel construction.
- ◆ Continuous Ice Flakes Output.
- ◆ CFC free compressor.
- ◆ Low water level Indication.
- ◆ Overload protection.
- ◆ Noiseless operation.
- ◆ Front & Back Ventilated Design.



HSN : 84186990

SPECIFICATIONS

PRODUCT NAME	ICE FLAKER		
MODEL	IF30D	IF50D	IF100D
Capacity of Ice flake (kg per 24 hrs.)	30 kg	50 kg	100 kg
Storage Capacity (kg ~ app.)	7 kg	15 kg	25 kg
Display	RED LED		
Storage bin	ABS Plastic Bin		
Ambient temperature	Below 25 °C for optimum working		
Construction	Stainless Steel		
Noise level	≤52 dB		
Water Supply	3 / 4 inch BSP		
Cooling Method & Refrigerant	Air & Refrigerant R290A		
Power supply	AC 220V / 50Hz		
Input Power	180W	280W	420W
Packing Dimension LxWxH	300X493X547mm	380X543X72mm	548X611X883mm
Gross Weight	36kg	45kg	64kg

* Specifications are subject to change without Notice

Monotech Lab Instruments and Calibration

B-9, Shabnam Building, MEK Industrial Estate, Opp. Fire Brigade, Mumbra, Dist. Thane 400612 Maharashtra