

FEATURES

MIM15

- ◆ Housing Material Stainless Steel Type
- ◆ R600A Refrigerant .
- ◆ Water Inlet Manually .
- ◆ Bullet Shape Ice.
- ◆ Low Water Alarm.
- ◆ Ideal Working Temperature 20°C.
- ◆ Transparent Window on Top Cover to view inside.



HSN : 84186990

SPECIFICATIONS

PRODUCT NAME	ICE MAKER
MODEL	MIM15
Ice capacity	15kg / 24h
Function	Ice Making
Storage Bin	2kg
Refrigerant	R600A
Ice Shape	Bullet shape
Housing material	Stainless steel type
Water Inlet method	Manually
Ice output each cycle	14 pieces
Ice Making Time	10 - 12 minutes in Ideal Conditions
Input power	150W
Power supply	AC 220V/50Hz
Packing Dimension (L x W x H)	375 x 320 x 465mm
Net / Gross Weight	12.5Kg / 14.5kg

* Specifications are subject to change without Notice

Monotech Lab Instruments and Calibration

B-9, Shabnam Building, MEK Industrial Estate, Opp. Fire Brigade, Mumbra, Dist. Thane 400612 Maharashtra