

LABORATORY WEIGHTS

E1 Class, E2 Class, F1 Class & F2 Class

Standard Weight Box

Class	Denomination	Number of weights
E1	1 mg to 200 gm	23 Pcs
E2	1 mg to 200 gm	23 Pcs
F1	1 mg to 200 gm	23 Pcs
F2	1 mg to 200 gm	23 Pcs
M1	1 mg to 200 gm	12 Pcs

Standard Weights

DENOMINATION	NUMBER OF WEIGHTS
500 GM	1 PC
1 KG	1 PC
2 KG	1 PC
5 KG	1 PC
10 KG	1 PC
20 KG	1 PC
50 KG	1 PC

Specifications :-

- Material : Stainless Steel 304 (Non-Magnetic)
- Density : 7.90 cm³
- Magnetic Susceptibility : <0.005
- Shape : Cylindrical with Knob

***All Weighing according to OIML - R111 Standards**

We can provide Calibration Certificate from NABL Accredited Lab (Charges applicable).

Monotech Lab Instruments and Calibration

B-9, Shabnam Building, MEK Industrial Estate, Opp. Fire Brigade, Mumbra, Dist. Thane 400612 Maharashtra India.

☎ +91-7738246010 / 8433872661

✉ info@monotechlabinstruments.com

🌐 www.monotechlabinstruments.com

