

FEATURES

MOS100

- ◆ Digital LCD Display.
- ◆ 100 ~ 2500rpm speed.
- ◆ Maximum Volume 20L.
- ◆ Stirring has three modes - Clockwise, Anticlock wise & Alternating.
- ◆ Stainless Steel Head.
- ◆ Digital type bushless motor.



HSN : 84798200

SPECIFICATIONS

PRODUCT NAME	OVERHEAD STIRRER
MODEL	MOS100
Max. Volume	20L
Display	LCD
Speed(r.p.m)	100~2500
Working Time	Continuous
Structure Exterior	Plastic shell
Motor power	60W
Setting & Speed mode	Digital
Timer	Yes
Max Viscosity	10000 Mpas
Plate size	200 × 300 mm
Plate load	50kg
Power Supply	Using Adaptor, AC 220V / 50-60Hz
Dimension & Weight	400 × 250 × 210 mm & 7.5kg

* Specifications are subject to change without Notice

Monotech Lab Instruments and Calibration

B-9, Shabnam Building, MEK Industrial Estate, Opp. Fire Brigade, Mumbra, Dist. Thane 400612 Maharashtra