

FEATURES

MUC3, MUC6, MUC9, MUC12, MUC25

- ◆ Inner Tank & Outer body made of Stainless Steel.
- ◆ Digital Control.
- ◆ RED LED display.
- ◆ Heating Temperature RT+5 ~ 80°C.
- ◆ Temperature with timer function.
- ◆ Drainage and cool fan available in 6L to 24L.
- ◆ LID & Basket as Standard Supply.



SPECIFICATIONS

HSN : 84562000

MODEL	MUC3	MUC6	MUC9	MUC12	MUC25
Tank Capacity (Approx.)	2.5L	5L	8L	11L	24L
Tank Size (mm)	240x140x100	300x155x150	300x240x150	330x300x150	500x300x200
Overall Size (mm)	270x170x240	330x180x310	330x270x310	360x330x310	550x330x360
Ultrasonic Power	100W	150W	200W	300W	500W
Heating power	100W	300W	300W	400W	500W
Tank Drain	NA	Drain Valve			
Time Setting	1-99 min.				
Heating Temp.	RT. +5 ~ 80°C				
Display	Digital RED LED				
Frequency	40kHz				
Lid Material & Inner Tank Material	SS 304				
Outer Tank Material	SS 201				
Power Supply	AC 220V / 50 Hz				
Weight (approx.)	4 kg	7.5 kg	9 kg	12 kg	17 kg

Monotech Lab Instruments and Calibration

B-9, Shabnam Building, MEK Industrial Estate, Opp. Fire Brigade, Mumbra, Dist. Thane 400612 Maharashtra