

FEATURES

MDV62T

- ◆ Touch screen Display.
- ◆ Memory upto 30 sets.
- ◆ Temperature Monitoring.
- ◆ Automatic Detection.
- ◆ Selectable Speed.
- ◆ Date & Time.
- ◆ Product Name Input.
- ◆ Optional Accessory - Low Viscosity Adaptor (0 Spindle).



HSN : 90278910

SPECIFICATIONS

PRODUCT NAME	DIGITAL VISCOMETER
MODEL	DV62T
Display	Touch screen
Measuring range	50-20,00,000 mPas
Measuring Accuracy	±2% (Newtonian liquid)
No of speed	8 Speed Option (0.3, 0.6, 1.5, 3, 6, 12, 30, 60 rpm)
Ambient Temperature	10°C to 50°C
No of Spindle(rotor)	4
Zero (0) Spindle(rotor)	Optional 0 Spindle (or 1 to 50mPas)
Interface	RS232C
Memory	30 sets
Date & Time	Yes
Single Measurement Auto Stop	Afer measurement result is stable, spindle stop automatically
Power supply	AC 220V/50Hz
Dimension (LxWxH) & Weight	475 × 200 × 375 mm & 6.8 kg

* Specifications are subject to change without Notice

Monotech Lab Instruments and Calibration

B-9, Shabnam Building, MEK Industrial Estate, Opp. Fire Brigade, Mumbra, Dist. Thane 400612 Maharashtra

VISCOMETER - TOUCH SCREEN VIEW

MDV62T

