

- ◆ Operating Mode - 2 types - Auto run sensor, Continuous run.
- ◆ Control through Knob.
- ◆ Bushless Motor.
- ◆ Supply with Standard Head.
- ◆ Soft silicone foot pads.
- ◆ Super shockproof.
- ◆ Aluminium Alloy Base.
- ◆ Stable body.
- ◆ Head for Vortex Mixer (**Optional**)
- ◆ Size 4 x 50ml, 8 x 10/15ml, 12 x 1.2ml/2m  
8 x 5ml, 24 x 0.5ml, 48 x 0.5ml
- ◆ Universal Plate , Flat gasket



**SPECIFICATIONS**

HSN : 84798200

PRODUCT NAME	VORTEX MIXER
MODEL	MVM60
Speed	100-2800rpm
Motor	Brushless DC motor
Adjustment method	Knob
Operation mode	2 Type - Auto run sensor, continuous run
Base material	Aluminium alloy
Shaking Movement	Orbital
Orbital diameter	4mm
Display	NA
Working Temperature	5°C-40°C
Power	60W
Power Supply	Using Adaptor, AC 220V/50Hz
Packing Dimension (L x W x H)	250x220x190mm
Net / Gross weight	3kg / 3.5kg

\* Specifications e subject to change without Notice